


Meiji Techno Boom Mounted Optics Part Numbers & Specifications



System# sunimac

Short UNIMAC are short-back focus, direct-to-camera, zoom optics. Small, easy to place and maneuver, Short UNIMAC is perfect for Quality Control, Inspection, Machine Vision and hundreds of other applications.

Part Numbers & Specifications

Part Numbers		Description
A	S-4300	Universal Boom Stand, requires FS focus block, vertical pole: 400mm, horizontal: 535mm. Base: 265mm x 265mm x 25mm.
B	No Eyepieces	No Eyepieces for this configuration
C	MS-50	Short UNIMAC, 0.7x to 4.5x zoom with Cmount Adapter included
D	MS-52/05	Inclinable coarse focus block with 5/8" pin, accepts S-4300 and S-4400 boom stands
E	FT190/115	150 watt fiber optic light source and power supply
F	CC2100	USB2 Digital Video Camera, 1.4MP, complete with software



**MEIJI
TECHNO**

Meiji Techno America 3010 Olcott Street Santa Clara, CA 95054
Toll Free USA: 1(800) 832-0060 (408) 970-4799 FAX: (408) 970-5054
Visit us online at www.meijitechno.com