



# ML9700 Series

## POLARIZING MICROSCOPES

### Features and Specifications

**Viewing head:** Binocular or trinocular body, 30° inclination, by model. Attachment fitting slotted for correct orientation.

**Revolving Nosepiece:** Ball-bearing quadruple nosepiece with one fixed and three centerable objective apertures.

**Eyepieces:** HWF10X Focusing eyepiece with cross-line reticule in right-hand fixed eyetube and a standard HWF10X Eyepiece in left-hand adjustable eyetube.

**Objectives:** Strain free objectives SM Plan 4X N.A. 0.10, SM Plan 10X N.A. 0.25, SM plan 40X N.A. 0.65

**Analyzer:** Rotatable analyzer in slider mount with 360° graduations

**Bertrand lens:** Bertrand lens with field-limiting stop, sliding for E-W insertion.

**Polarizer:** Swing-out click/clamp mount, fully rotatable with click stop at 0° and 90°

**Condenser:** Strain free Achromatic condenser, N.A. 1.25 with aperture iris diaphragm and swing-out filter holder.

**Illuminator:** Built-in base Koehler illuminator, 50 watt halogen lamp with intensity control.


### Models:

**ML9720** Binocular Polarizing Microscope

Size and Weight 200mm(W) x 270mm(D) x 435mm(H)  
8.6Kg

**ML9730** Trinocular Polarizing Microscope

Size and Weight 200mm(W) x 270mm(D) x 480mm(H)  
8.8Kg



ML9720

Agent :



**MEIJI TECHNO CO.,LTD.**

322-1, Chikumazawa Miyoshi-machi, Iruma-gun  
Saitama 354-0043, Japan

E-mail : [meiji@meijitechno.co.jp](mailto:meiji@meijitechno.co.jp)

Phone : (0)49-259-0111 Fax : (0)49-259-0113

**MEIJI TECHNO AMERICA**

2186 Bering Drive, San Jose, CA., 95131, USA

E-mail : [sales@meijitechno.com](mailto:sales@meijitechno.com)

Phone : 408-428-9654 Fax : 408-428-0472

Toll free : 800-832-0060

<http://www.meijitechno.com>

